Quick Guide

AgraStrip[®] Pro WATEX[®] AgraVision™ Pro

The AgraStrip[®] Pro WATEX[®] test system enables the rapid and simple on-site quantification of mycotoxins in a variety of agricultural commodities. Both extraction and assay have been streamlined, reducing steps to a bare minimum while setting new standards in accuracy and usability.

AgraStrip® Pro Watex® mycotoxin test kits include:

- 40 test strips
- 40 extraction buffer bags
- 40 cartridges for AgraVision™ Pro Reader
- 1 bottle of 50 ml assay dilution buffer
- 40 yellow or white pipette tips

- 40 blue pipette tips (80 for T-2/HT-2)
- 40 microcentrifuge tubes (dilution tubes) (80 for T-2/HT-2)
- 40 Whirl-Pak® filter bags
- 1 tube holder
- \cdot 1 QR code card

Item	Description	ltem No.
AgraStrip® Pro Total Aflatoxin WATEX®	40 tests	10006322
AgraStrip® Pro Total Fumonisin WATEX®	40 tests	10006318
AgraStrip® Pro Deoxynivalenol WATEX®	40 tests	10006320
AgraStrip® Pro Ochratoxin A WATEX®	40 tests	10006332
AgraStrip® Pro Zearalenone WATEX®	40 tests	10006316*
AgraStrip® Pro T-2/HT-2 WATEX®	40 tests	10006614
AgraVision™ Pro	Reader	10006324
AgraVision™ Pro Printer	Printer	10006483
AgraVision™ Pro Scanner	Scanner	10006484
AgraVision™ Pro Wi-Fi dongle	Wi-Fi dongle	10006623

* To perform test according to FGIS protocol, please use the kits which include material for 50 g sample extraction.







AgraStrip[®] Pro WATEX[®] – AgraVision[™] Pro

To ensure accurate results, please read the product manual before performing the assay.

Procedure on corn



